

NOAH 4.2.0

Patch Build 1750

May 2012

This patch can be installed on the following systems:

- NOAH 4.2.0 (build 1749)
- NOAH 4.2.0 (build 1747)

The patch installation will not run on any other Noah installations.

This patch addresses the following:

This patch will update Noah client and Noah System Server components. The following issues are addressed with this patch:

- Problem saving actions from a module. In some cases when a module saves an action Noah will fail. This happens in certain situations if the patient data has been updated prior to saving a new action.
Installing this patch ensure that Noah client handles the action and session list correct also after editing patient data.
- When a client PC returns from a suspend or sleep state it has been observed that actions are saved on another patient than intended.
Installing this patch ensure that Noah System database interface is updated and saves actions on the correct patient.
- When a client PC returns from a suspend or sleep state it has been observed that the client in some situations loses connection to the Noah server after a little while.

To ensure that the patch is installed properly:

The patch should be applied to both Noah Server and Noah client PC's. You can see in the Windows Add/Remove programs list that the patch has been installed. You will need to view installed upgrades.

After the patch has been installed the Noah System about box will list the Noah version: 4.2.0 (build 1750).